

TRAY RE-ORDERING MACHINE



The MT-S/D can handle these kind of containers:













MT is a single or double unit to collect vials into trays.

Quincux system allows vials collection into orderly rows which make trays ready for their transfer to a freeze-drier for **lyophilisation process**.

When the tray is full, the machine stops automatically. If it is a double tray loader when the first tray is full the vials are directed to the second tray and vice versa. This enables continuous operation.

TECHNICAL DATA	MT-S	MT-D	
DEVELOPMENT	Linear	Linear	
MOVEMENT	Intermittent	Intermittent	
OUTPUT	Up to 150 pcs/min	Up to 300 pcs/min	

