

# MDT SERIES

DEPYROGENATION TUNNELS FOR AMPOULES AND VIALS  
FOR SMALL VOLUME PARENTERAL (SVP)



- Hot chamber with dedicated make-up unit and air-flow control.
- Automatic balancing and air flow adjustment.
- Pressure cascade according to process requirements.
- Outfeed chamber with dedicated cooling coil and 100% air recycling.
- Dry-Heat sterilization of outfeed chamber.

#### TECHINCAL DATA

<b>PRODUCT RANGE</b>	<b>Vials</b>	2 R / 500 ml
	<b>Cartridges</b>	1 ml / 3 ml
	<b>Ampoules</b>	1 ml / 30 ml
<b>OUTPUT</b>	<b>Pieces/m</b>	According tunnel size and containers to be processed