

MW SERIES

ROTARY CONTINUOUS MOTION MACHINES FOR INTERNAL WASHING AND DRYING OF AMPOULES AND VIALS FOR SMALL VOLUME PARENTERAL (SVP).



- 8 internal and 3 external washing/blowing stations
- Pre-washing in Ultrasonic Bath
- Water Recycling
- Automatic drainage and drying
- Container siliconizing

TECHINCAL DATA

PRODUCT RANGE	Vials	2 R / 500 ml
	Cartridges	1 ml / 3 ml
	Ampoules	1 ml / 30 ml
OUTPUT	Pieces/m	50 / 400